

Portion cupboard Metos Cube-4

Cupboard to accommodate patients meal trays and hold them at right temperatures for extended time after normal meal distribution. To be used as complement to any Burlodge meal tray delivery system

- Ideal in areas with flexible or variable serving times
- Maintains meal components at right hot and cold serving temperatures
- Avoids wasted meals due to storage in unsafe conditions
- Minimize late trays and trolley returns to the kitchen

Features:

- isolated thermal barrier divides the flat tray in a hot and a cold zone
- digital display, automatic start timer
- tray: Metos Burlodge 530 x 325 mm or 575 x 325 mm
- small fot print: external measures (L \times D \times H): 681 \times 620 \times 475 mm (681 \times 620 \times 1222 with stand)
- material: brushed stainless steel, recyclable insulation, tempered glass doors
- several units stackable on top of each other to add capacity

Accessories

- stand 4216112
- stacking spacers 4216111





Portion cupboard Metos Cube-4

Product capacity	4 trays, distance 80 mm
Item width mm	681
Item depth mm	620
Item height mm	475
Package volume	2.275
Unit of volume	hl
Package volume	2.275 hl
Package length	70
Package width	65
Package height	50
Package unit of dimension	cm
Package dimensions (LxWxH)	70x65x50 cm
Net weight	40
Net weight	40 kg
Gross weight	50
Package weight	50 kg
Unit of weight	kg
Connection power kW	0.6
Fuse Size A	3
Connection voltage V	230
Number of phases	1N
Frequency Hz	50
Type of electrical connection	Plug